

Historic, Archive Document

Do not assume content reflects current scientific knowledge, policies, or practices.

Reserve
aQC452
.5
.R69
1977

FINAL REPORT

INSTRUMENTATION FOR MEASURING IRRADIANCE IN THE UV-B REGION

J. D. Rowan
K. H. Norris

Instrumentation Research Laboratory
Agricultural Marketing Research Institute
Beltsville Agricultural Research Center
Beltsville, Maryland 20705

EPA-IAG-D6-0168

Project Officer:

R. J. McCracken
Agricultural Research, Science and Education Administration
U.S. Department of Agriculture
Washington, D.C. 20250

Prepared for
Environmental Protection Agency
BACER Program
Washington, D.C. 20460

**United States
Department of
Agriculture**



National Agricultural Library

Annual Report 1977

EPA Interagency Program on Biological and Climatic Effects (BACER) Instrumentation for Measuring Irradiance in the UV-B Region

SUMMARY

The responsibility of the Instrumentation Research Laboratory (IRL) in this program was developing portable instruments for use by biologists to measure UV-B irradiance in growth chambers, greenhouses, and field plots. A simple UV-B radiometer and two UV-B spectroradiometers have been designed, constructed, tested, and put into use in the UV-B research program. Each of these instruments is now being manufactured by commercial firms.


The two spectroradiometers differ only in the monochromators: one has a single holographic grating; the other, two holographic gratings for greater stray light rejection. The spectroradiometers automatically scan the 250-400 nm region in less than 5 minutes, printing a tape of the corrected irradiance as a function of wavelength. The input is cosine corrected by a specially designed teflon-bubble diffuser coupled to the input slits. The output of the monochromator is measured with a multiplier-type phototube and a logarithmic response amplifier. The amplifier output is digitized with a digital voltmeter, and the digital output is interfaced with a desk-top programmable calculator.

The design provided for the spectroradiometer system response to be stored in the calculator so that, as the spectrum is scanned, the calculator corrects the measured signal for instrument calibration and outputs the true spectral irradiance of the source being measured.

The programmable calculator controls the operation of the spectroradiometer so that, on command, scanning is initiated and readings are recorded for each nanometer interval. The calculator prints the wavelength and irradiance for each wavelength interval and, at the end of the scan, reverses the wavelength drive and returns the monochromator to the starting wavelength. The calculator is programmed to print an integral for a programmed action spectrum at the end of each scan. At the completion of the scan, the data can be stored on a magnetic tape, if desired, for future analyses.

Provision is included for a precise check of wavelength: the spectrum of a miniature low-pressure mercury-arc lamp is scanned, and the calculator computes the position of the 253.7 nm and the 296.7 nm lines to a precision of ± 0.01 nm.

The performance specifications approach the requirements for UV-B measurements stated in NBS publication #20, "Optical Radiation News," dated April 1977.



DEPT OF AGRICULTURE
JAN 21 1936

Annual Report 1977

EPA Interagency Program on Biological and Climatic Effects (BACER) Instrumentation for Measuring Irradiance in the UV-B Region

Our UV-B measuring instruments developed for this project meet the requested requirements. The broadband radiometer is small, hand-held, battery operated, and has fast response for rapid measurements. The spectroradiometers have the following features:

- (1) Fast, accurate, and reproducible.
- (2) Convenient, automatic single-key-stroke operation.
- (3) Programmable calculator-controlled scanning, Fig. 1.
- (4) Logarithmic amplifier for wide dynamic range.
- (5) Calculates and prints true spectral irradiance for each nm wavelength, Fig. 2.
- (6) Calculates and prints weighting function ($A_{\Sigma 9}$).
- (7) Programmed-wavelength-calibration check.
- (8) Double monochromator unit has very low stray light.
- (9) Bandwidth of 2 nm.
- (10) Temperature stable.
- (11) Receptor has excellent cosine response.
- (12) Spectra are stored on magnetic tape for efficient data processing.
- (13) Portable--operates on small lab cart.
- (14) Minor disturbance of area by the "measuring head."

Our spectroradiometer performance specifications approach the requirements for UV-B measurements stated in NBS publication #20, "Optical Radiation News," dated April 1977. These requirements are cited in attached copy of "Making a UV Measurement?," Electro-Optical Systems Design, 9(6):17, 1977. (page 7)

Cosine Response for UV-B

The spectroradiometers, as well as other UV-B radiation meters developed in our program, incorporate a new design for cosine-corrected input optics. We tested all available cosine-correction schemes, including integrating spheres, diffusing reflectors, and sintered-quartz diffusers. None of these gave adequate performance for the 250 to 370 nm region, so we developed the

teflon bubble diffuser. The spectral transmission of this teflon diffuser is shown in Fig. 3. Fig. 4 is a drawing of the input optics for the spectroradiometers. We supplied a similar receptor to our cooperators that have a Gamma Scientific monochromator, Fig. 5. The response of this type of diffuser to radiation from different angles is compared with that of a commercial instrument in Fig. 6. The teflon bubble receptor provides excellent cosine correction, and has a stable surface which can be readily cleaned. This type of diffuser is now used in commercial instruments.

Wavelength Accuracy

To check spectroradiometer wavelength accuracy, a low-pressure mercury-arc lamp is scanned and the mercury line location computed. The wavelengths of the 253.7 and the 296.7 nm lines are measured with a readout precision of ± 0.01 nm, so that wavelength shifts as low as 0.02 nm are readily detected, Fig. 7. This field check can be made quickly and routinely between measurements if desired.

Spectroradiometers

Our UV-B spectroradiometers are now used routinely by technicians to measure UV-B irradiance in growth chambers, greenhouses, and field plots. These instruments are identified by the acronyms: IRLSpec-S for the single monochromator version; IRLSpec-D for the double monochromator; and IRLSpec-SO for the commercial model single monochromator. Numbers shown on the graphs--for example, 1770815.09--identify instrument (first digit), year (next two digits), month, day, and scan number for that day. The first digit may be 1, 2, or 3 for IRLSpec-S, IRLSpec-D, or IRLSpec-SO, respectively.

Measured performance specifications for the IRLSpec-S spectroradiometer are as follows:

- Wavelength range - 250 to 370 nm
- Readout interval - 1 nm
- Scanning speed - 0.5 nm/sec.
- Spectral bandpass - 2 nm
- Wavelength reproducibility - ± 0.02 nm
- Wavelength accuracy - ± 0.05 nm at 296.7 nm
- Radiometric reproducibility - $\pm 2\%$
- Radiometric accuracy - $\pm 5\%$
- Radiometric range - 0.001 to $2000 \text{ mW} \cdot \text{m}^{-2} \cdot \text{nm}^{-1}$
- Stray light - less than 2×10^{-4} at 285 nm as tested
with a xenon arc lamp filtered with 0.5-mm-thick
cellulose acetate
- Cosine correction - within $\pm 10\%$
- Output - printed tape direct reading in wavelength
and corrected irradiance at each nanometer interval
from 250 to 370 nm
- Size of measuring head - less than $28 \times 20 \times 10$ cm
- Temperature stability less than 0.04 nm wavelength shift
for 25°C temperature change.

The IRLSpec-D spectroradiometer, with double-grating monochromator, provides better stray light rejection. This unit incorporates all the features of the single monochromator unit although the measuring head is, of necessity, slightly larger. Observed wavelength instability of the double monochromator was caused by grating-sync-drive cable temperature sensitivity. Our technician, George Button, solved this temperature instability with his cable-stringing technique. The grating-drive cable spring was eliminated, and the cables were crossed to provide automatic temperature compensation and wavelength stability.

Measured performance specifications for the IRLSpec-D spectroradiometer are as follows:

- Wavelength range - 250 to 370 nm
- Readout interval - 1 nm
- Scanning speed - 0.7 nm/sec.
- Spectral bandpass - 2 nm
- Wavelength reproducibility - ± 0.02 nm
- Wavelength accuracy - ± 0.05 nm at 296.7 nm
- Radiometric reproducibility - $\pm 2\%$
- Radiometric accuracy - $\pm 5\%$
- Radiometric range 0.0001 to 1000 $\text{mW}\cdot\text{m}^{-2}\cdot\text{nm}^{-1}$
- Stray light - less than 5×10^{-8} at 285 nm as tested
with a xenon arc lamp filtered with 0.5-mm-thick
cellulose acetate
- Cosine correction - within $\pm 10\%$
- Output - printed tape direct reading in wavelength and
corrected irradiance at each nanometer interval from
250 to 370 nm
- Size of measuring head - less than 28 x 25 x 10 cm
- Temperature stability - less than 0.06 nm shift for
25°C temperature change.

The IRLSpec-D spectroradiometer was field tested at the University of Colorado's Snowmass, Colorado, high-altitude site during August. The instrument was disassembled into components and securely packed into a regular hard-side luggage as a shipping container and transported as "checked" baggage to and from Colorado. Spectral measurements were made at three elevations, 2377 m (7800 ft.), 2980 m (9777 ft.), and 3452 m (11326 ft.). This necessitated disassembling and repacking the instrument and transporting it up and down the mountain on unimproved roads to each site. Upon return to Beltsville, the instrument continued to perform with only minor adjustments.

The commercial version of our IRLSpec-S, single monochromator spectroradiometer was manufactured to our specifications. This IRLSpec-SO required constant operator assistance, but our circuit modifications converted it to automatic operation with single-key stroke to scan the spectrum, return, and reset.

Measured performance specifications for this commercial IRLSpec-SO spectroradiometer are as follows:

Wavelength range - 250 to 370 nm
Readout interval - 1 nm
Scanning speed - 1.0 nm/sec.
Spectral bandpass - 2 nm
Wavelength reproducibility - ± 0.02 nm
Wavelength accuracy - ± 0.05 nm at 296.7 nm
Radiometric reproducibility - $\pm 2\%$
Radiometric accuracy - $\pm 5\%$
Radiometric range - 0.001 to 200 $\text{mW}\cdot\text{m}^{-2}\cdot\text{nm}^{-1}$
Stray light - less than 1×10^{-4} at 285 nm as tested
with a xenon arc lamp filtered with 0.5-mm-thick
cellulose acetate
Cosine correction - within $\pm 10\%$
Output - printed tape direct reading in wavelength and
corrected irradiance at each nanometer interval from
250 to 370 nm
Size of measuring head - less than 25 x 21 x 13 cm
Temperature stability - less than 0.04 nm shift for
25°C temperature change.

Spectroradiometer Sensitivity and Dark Current

The photomultiplier-tube dark current, system correction, and standard-lamp spectral irradiance are plotted in Fig. 8 for the IRLSpec-S with Hamamatsu R166 phototube (solar blind response) and Corion SF-1.25 solar blind filter. The IRLSpec-D photomultiplier dark current (Hamamatsu R212, S-5 response and Corion SF-1.25 solar blind filter), system correction, and standard lamp spectral irradiance are plotted in Fig. 9. Plots of the IRLSpec-SO photomultiplier-tube dark current (Hamamatsu R166, solar blind response, and Corning 7-54 filter), system correction, and standard lamp spectral irradiance are shown in Fig. 10.

The instruments developed in our program are now manufactured. A broadband UV-B radiometer is available from Optronics Labs, Inc. Spectroradiometers similar to ours are advertised in Optical Spectra, 12(1):13,51, and 60, dated January 1978, and are available from three firms:

EG&G, Inc.
35 Congress Street
Salem, Massachusetts 01970

Gamma Scientific, Inc.
3777 Ruffin Road
San Diego, California 92123

Optronics Labs, Inc.
7676 Fenton Street
Silver Spring, Maryland 20910

Irradiance Spectra

Typical spectra from our IRLSpec-D are plotted in Fig. 11, 12, 13, and 14; from our IRLSpec-S, Fig. 15 and 16; and from our IRLSpec-SO, Fig. 17.

Irradiances of FS40 and FBZS40 lamps are compared in Fig. 11. A typical curve for an FS40 lamp with and without a cellulose acetate (CA) filter is shown in Fig. 12. CA filters with FS40 lamps are compared in Fig. 13. Measured irradiances of four FS40 lamps used to enhance the sun's UV-B at Snowmass, Colorado's 2980 m (9777 ft.) site, are plotted in Fig. 14. The fixtures were two lamps mounted in pairs at a Z of 110 cm and filtered with 0.127-mm (0.005-in.) CA. The change in irradiance of a FBZS20 lamp after aging 16.5 hr. is shown in Fig. 15. The irradiance measurements of two Sylvania F15T8 CW lamps at Z = 8 cm are compared for the IRLSpec-S, Fig. 16; the commercial IRLSpec-SO, Fig. 17; and the IRLSpec-D, Fig. 18. Irradiance for a 15-watt daylight lamp is shown in Fig. 19.

Spectral data of the sun at Beltsville, Maryland, during the year are plotted in Fig. 20 and 21. The sun curve on a very clear day illustrates the stray light at wavelengths below 289 nm for the IRLSpec-D, Fig. 22. The action spectra of measured Beltsville sunshine with weighting function 9 are shown in the region of 294 to 318 nm at the bottom of Fig. 20.

The weighting function plotted on the graphs is the A₂₉ equation developed in cooperative research at BARC:

$$\hat{y} = \left[0.25(\lambda/228.178)^{9.0} \right]^4 \exp \left[4 - (\lambda/228.178)^{9.0} \right]$$

The equation was used to compute UV-B sun equivalent (UVSEB), indicated on the graphs. Details of its development are presented in the Agricultural Equipment Laboratory report.

Spectral irradiance data are now plotted automatically with the HP9815A calculator interfaced directly and controlling a HP9872A plotter. After the data are recorded on the cassette, selected scans are loaded into the calculator; and the plotter draws and labels the graph to programmed dimensions, and then automatically plots the data (Fig. 11, 12, 13, and 14). A curve is plotted from 250 nm to 370 nm in 25 seconds.

Narrow-band Radiometers

Two narrow-band, probe-type radiometers were commercially manufactured to our specifications. These radiometers have solar-blind photo-multiplier-type detectors (Mfg. No. R166), a peak response at 294 nm and 299 nm, respectively, and a half-band width of 5 nm. Irradiance readout is a 3 1/2 digit display located in the remote electronic unit with switched ranges from 10⁻⁴ to 10⁻¹¹ watt·cm⁻²·nm⁻¹. The small sensor head (5x5x16 cm), with teflon dome receptor, is attached by 2.75-m cable for easy placement within growth chambers.

Broad-band UV-B Radiometers for Routine Monitoring

A small battery-powered radiometer, the IRLMeter, was designed and constructed for measurement of total UV-B radiation from artificial sources. The instrument has a solar-blind vacuum photo-diode (Hamamatsu R403), an integrated circuit amplifier, and a microammeter packaged in a meter case. The circuit, as shown in Fig. 23, provides for four decades of range switching. The relative wavelength response, as shown in Fig. 24, peaks at 300 nm and is relatively flat from 280 to 320 nm and the radiometer is not sensitive to radiation longer than 400 nm. Typical comparative data of the IRLMeter with the IRLSpec-S spectro-radiometer are shown in Fig. 25.

This broad-band UV-B radiometer sensitivity is adjusted to a full scale reading of 2.0 UV-B sun equivalents, Beltsville (UVSEB) when measuring the output of FS40 fluorescent lamps filtered with 0.127-mm (0.005 in.) CA. Correction factors were developed for use with FS40 lamps with 0.254-mm (0.010 in.) CA filter, and for use with BZS fluorescent lamps.

The IRLSun-meter radiometer is a modification of the IRLMeter. The new calibration reads directly in sun equivalents, Beltsville (UVSEB AΣ9). Calibrations for two suns (UVSEB AΣ9) full scale were developed for four sources:

- (1) Sunshine for a clear day (6/8/77), 1:00 p.m. EDT
- (2) FS40 lamps filtered with 0.127-mm (0.005 in.) CA
- (3) FS40 lamps filtered with 0.254-mm (0.010 in.) CA
- (4) FBZS40 CLG lamps (unfiltered)

A commercial version of the IRLMeter has been manufactured, the Optronic Model 725. These Model 725 radiometers were calibrated by IRL for our cooperating laboratories from IRLSpec-D irradiance values for five sun equivalents full scale, with FS40 source filtered by 0.127-mm (0.005 in.) CA aged 6 hr. A calibration scale factor was developed for 0.191-mm (0.0075 in.) CA, 0.254-mm (0.010 in.) CA, 1.52-mm (0.036 in.) petri dish, and FBZS40 and FBZS20 lamps.

An Optronic Model 725 broad-band radiometer was calibrated from IRLSpec-D irradiance values for 10 UVSEB full scale with a Rayonet F8T5 RPR3000A (8 watt) lamp source filtered with 0.127-mm (0.005 in.) CA supplied by the Peoria cooperators and then aged 6 hr. This Peoria CA filter absorbed less at shorter wavelength than Beltsville CA, Fig. 26, 27, and 28, significantly changing the AΣ9 action integral.

Rayonet F8T5 RPR3000A (8 watt) lamps have strong energy in the 254 nm region, the energy approaching the peak value at 313 nm (ratio 1.4). The peak energy of the FS40 at 313 nm is approximately 700 times the energy at 254 nm. This 254 nm energy can be removed with CA filter, but it probably accelerates the aging of the CA filter. The irradiance of the Rayonet RPR3000A with the CA filter at 20 cm approximates the energy of the FS40 with CA filter at 50 cm and 75.1 cm, Fig. 29, 30, 31, and 32.

Eye Protectors

Personal eye protectors were evaluated for their UV-B attenuation, and two FS40 lamps at 50 cm from the IRLSpec-D spectroradiometer receptor were used. An aperture through black cloth was necessary to prevent significant leakage between lens and receptor, even though they were placed as close as possible. Plots of plastic goggles No. 6, 7, and 8 (dashed line) are at the 0.001 irradiance level in Fig. 33. No. 5 was a pair of plastic goggles, and No. 2, a plastic face shield. UV-B blocking of prescription eyeglasses (glass and plastic) and sunglasses are shown in Fig. 34.

EOSD

Electro-Optical Systems Design

The Engineering Magazine of Electro-Optical and Laser Technology
June 1977, Vol. 9, No. 6

MILTON S. KIVER
Publisher

STEPHEN SKILNYK
Co-Publisher & Vice President

RICHARD CUNNINGHAM
Editor

WILLIAM S. HUDSPETH
Associate Editor

WILLIAM D. ASHMAN

LOIS ESMAIL
Circulation Manager

GENE BRIESKE
Marketing Services

EDITORIAL ADVISORY BOARD

L. Beiser, *Consulting Physicist*
D. Belforte, *Avco Everett Research Lab*
R. Buzzard, *Lexel*
D. Casasent, *Carnegie-Mellon Univ.*
W. Hunter, *Naval Research Lab*
G. Klauminzer, *Molelectron*
H. Lavin, *General Electric*
D. Lockie, *Watkins-Johnson*
M. Pasturel, *Coherent Radiation*
W. Ruderman, *Interactive Radiation*
R. Watson, *EG&G*

EOSD is the sponsor of the
ELECTRO-OPTICS/LASER Conferences
and Expositions.

ELECTRO-OPTICAL SYSTEMS DESIGN.
Copyright © 1977 by Milton S. Kiver Publications, Inc. Published monthly by Milton S. Kiver Publications, Inc. Printed in U.S.A. Controlled circulation paid at Pontiac, Illinois and Chicago, Illinois 60607. Subscription price: \$3 per copy, \$20 per year. POSTMASTER: Send Form 3579 to ELECTRO-OPTICAL SYSTEMS DESIGN, 222 West Adams, Chicago, Ill. 60606. (312) 263-4866.



OFFICIAL PUBLICATION

JUNE, 1977

News

- 4 Laser "Writes" in Liquid Crystal
- 4 And Now: The Free Electron Laser
- 6 The Latest in Wrist Wear
- 6 Start Clicking for Nikon's Photo Contest
- 6 Heat-Generated Ultrasonic Waves
- 8 Future Brightens for LCDs
- 12 Quality Control Session to be Featured
- 12 Coherent Laser Radiation Shrinks to 38 nm
- 12 Unconventional Form for CRT Terminal
- 13 Need a Laser Diode?
- 14 Instant Movie Camera Debuts
- 17 BRH OKs First Laser Safety Standard Variance
- 17 Making a UV Measurement?

Making a UV Measurement?

Washington, D. C. — The increasing interest in ultraviolet optical measurements has spurred the National Bureau of Standards' Optical Radiation and Radiometric Physics groups to look into the matter. (A special chapter of the NBS Self-study Manual on Optical Radiation Measurements will be devoted to making state-of-the-art accuracy UV measurements and is planned for later this year.) Over the past few months, NBS has investigated and characterized several currently available spectroradiometers suited for use in the UV, in a program sponsored by the EPA.

To summarize their findings (see *Optical Radiation News*, April 1977), the current state-of-the-art for spectral irradiance measurements in the 250-350 nm range varies widely. The uncertainty for a highly controlled laboratory measurement of radiation of a simple character is estimated at about 3%. At the other extreme, the measurement of very complex radiation under very unfavorable field conditions could hit 26%. Typical

uncertainties for lab and field are estimated to be about 6% and 15%, respectively.

What can be done to improve these "best" UV measurements? NBS believes that what is needed most is a field-portable spectroradiometer with some impressive characteristics. This instrument would not only have negligible out-of-band leakage ($\approx 10^{-6}$) and high wavelength accuracy (≈ 0.03 nm) but would also be much less sensitive to temperature and mechanical shock, be hermetically sealed and have an easily replaceable window on its entrance port. Difficult field measurements could be made to within 10% uncertainty with such an instrument. NBS also suggests microprocessor control for this device.

(Ed. note: The technology for such a device exists. But such a spectroradiometer would necessarily be a "big ticket" item even if we disregard the mechanics of just how the development expenses are amortized. There is some question of whether the potential market is big enough to be able to absorb the development costs.) □

LIST OF FIGURES

- Fig. 1 - Spectroradiometer components
- Fig. 2 - Typical UV spectra
- Fig. 3 - Teflon spectral transmission
- Fig. 4 - IRLSpec-S receptor optics
- Fig. 5 - IRL receptor optics for gamma
- Fig. 6 - Receptor cosine response
- Fig. 7 - Hg line wavelength check - IRLSpec-D
- Fig. 8 - IRLSpec-S sensitivity and dark current
- Fig. 9 - IRLSpec-D sensitivity and dark current
- Fig. 10 - IRLSpec-SO sensitivity and dark current
- Fig. 11 - Lamp spectra (unfiltered) FS40, FBZS40CLG1076 and weighting function 9
- Fig. 12 - Lamp spectra FS40 and cellulose acetate 6.5 hr. Z = 30 cm, Z = 50 cm
- Fig. 13 - Lamp spectra FS40 and 0.127 mm, 0.191 mm, 0.254 mm cellulose acetate
- Fig. 14 - Lamp spectra, Snowmass, Colorado
- Fig. 15 - Lamp aging - BZS20CLG0377
- Fig. 16 - Sylvania F15T8 CW lamp spectra - IRLSpec-S
- Fig. 17 - Sylvania F15T8 CW lamp spectra - IRLSpec-SO
- Fig. 18 - Sylvania F15T8 CW lamp spectra - IRLSpec-D
- Fig. 19 - GE F15T8D lamp spectra - IRLSpec-D
- Fig. 20 - Sun spectra - Beltsville, MD
- Fig. 21 - Sun spectra - Beltsville, MD - IRLSpec-S
- Fig. 22 - Sun - Beltsville, MD - IRLSpec-D
- Fig. 23 - UV-B IRLMeter circuit
- Fig. 24 - IRLMeter spectral response
- Fig. 25 - Regression IRLSpec-S vs. IRLMeter

- Fig. 26 - Rayonet F8T5 RPR 3000A^o spectra new cellulose acetate, Z = 16.3 cm
- Fig. 27 - Rayonet F8T5 RPR 3000A^o spectra cellulose acetate aged 6 hr., Z = 16.3 cm
- Fig. 28 - Rayonet F8T5 RPR 3000A^o spectra cellulose acetate aged 6 hr., Z = 20 cm
- Fig. 29 - Rayonet F8T5 RPR 3000A^o and sun (Beltsville) spectra
- Fig. 30 - Spectra FS40 and Rayonet F8T5 RPR 3000A^o
- Fig. 31 - Spectra FS40 and cellulose acetate aged 0.5 hr., Z = 50 cm
- Fig. 32 - Spectra FS40 and cellulose acetate aged 6 hr., Z = 50 cm
- Fig. 33 - Spectra - eye protective device
- Fig. 34 - Spectra prescription eyeglasses and sunglasses





Fig. 1. Spectroradiometer components

FILE -17 **D 3**
 TYPE 2
 USED 1208
 MAX 1208
F540 Z=50 cm

2780203.09
 250000.0039
 251000.0048
 252000.0305
 253000.1084
 254000.2177
 255000.1071
 256000.0087
 257000.0030
 258000.0027
 259000.0027
 260000.0022
 261000.0024
 262000.0031
 263000.0033
 264000.0095
 265000.0389
 266000.0499
 267000.0438
 268000.0432
 269000.0652
 270000.1144
 271000.1838
 272000.2473
 273000.3516
 274000.5030
 275000.7620
 276001.0823
 277001.3153
 278001.6617
 279002.1661
 280003.0536
 281004.0457
 282004.7257
 283005.6125
 284006.8764
 285008.8472
 286011.0359
 287012.6380
 288014.3190
 289016.9874
 290020.5971
 291023.4209
 292025.7890
 293028.3414
 294031.4599

295035.6576
 296041.5418
 297044.7439
 298047.0842
 299049.3556
 300052.6190
 301057.2066
 302060.8413
 303063.8780
 304065.5058
 305067.2258
 306068.9599
 307071.0456
 308073.4446
 309074.8453
 310075.3275
 311078.3118
 312088.0010
 313095.7436
 314095.1361
 315078.5408
 316074.1841
 317073.4798
 318074.7080
 319073.8814
 320070.0362
 321067.2675
 322066.0442
 323066.6018
 324065.1532
 325060.8393
 326057.3499
 327055.6234
 328055.8803
 329054.1393
 330050.1182
 331046.6747
 332044.9649
 333045.0813
 334044.2513
 335040.8120
 336036.6656
 337034.9533
 338034.4372
 339032.7606
 340030.0634
 341027.7691
 342026.5048
 343025.9462
 344024.6910

345022.5070
 346020.5725
 347019.5569
 348019.0976
 349018.1677
 350016.6673
 351015.1844
 352014.3864
 353014.0250
 354013.3521
 355012.0807
 356011.0327
 357010.4609
 358010.2614
 359009.6726
 360008.8263
 361008.0685
 362007.6653
 363007.5352
 364014.7558
 365033.0470
 366023.8414
 367012.3363
 368007.0350
 369005.1119
 INT.9 FROM
 250
 80.2453
 SUNS = 26.2240
 SUM FROM
 250
 TO
 369
 MILLIWATTS/MSQ
 3428.0341
 SUM FROM
 270
 TO
 290
 MILLIWATTS/MSQ
 117.1260
 SUM FROM
 290
 TO
 320
 MILLIWATTS/MSQ
 1880.9138
 SUM FROM
 320
 TO
 340
 MILLIWATTS/MSQ
 1059.7177

Fig. 2. Typical UV spectra

0.020" TEFLON DISK




Fig. 3. Teflon spectral transmission




Fig. 4. IRLSpec-S RECEPTOR OPTICS

3-31-77





Fig. 5. IRL RECEPTOR OPTICS FOR GAMMA

3-31-77



11-3-76

Fig. 6. Receptor cosine response

H RUN SPECTRA
 B REVERSE
 D RECORD DATA
 E STOP
 F FORWARD
 G PRINT DATA
 I INT. 9
 J W. L. CHECK
 M SUMS
 N DARK ADJUST
 DATE
 2780203.00
 W. L. CAL.

2780203.01

2940002.164
 2950020.760
 2960043.339
 2970047.468
 2980029.692

MERCURY LINE AT
 296.738

W. L. CAL.

2780203.02

2940002.129
 2950020.865
 2960043.339
 2970047.444
 2980029.864

MERCURY LINE AT
 296.736

W. L. CAL.

2780203.03

2940002.155
 2950020.933
 2960043.419
 2970047.576
 2980030.002

MERCURY LINE AT
 296.738

W. L. CAL.

2780203.04

2940002.136
 2950020.808
 2960043.359
 2970047.510
 2980030.015

MERCURY LINE AT
 296.739

W. L. CAL.

2780203.05

2940002.098
 2950020.808
 2960043.319
 2970047.532
 2980030.210

MERCURY LINE AT
 296.742

W. L. CAL.

2780203.06

2940002.140
 2950020.665
 2960043.339
 2970047.488
 2980029.754

MERCURY LINE AT
 296.738

W. L. CAL.

2780203.07

2940002.192
 2950021.214
 2960043.419
 2970047.270
 2980029.117

MERCURY LINE AT
 296.731

Fig. 7. Hg line wavelength check - IRLSpec-D




Fig. 8. IRLSpec-S sensitivity and dark current




Fig. 9. IRLSpec-D sensitivity and dark current




Fig. 10. IRLSpec-60 sensitivity and dark current




Fig. 11 Lamp spectra (unfiltered) FS40, FBZS40CLG1076 and weighting function 9




Fig. 12 Lamp spectra FS40 and cellulose acetate 6.5 hr. Z = 30cm, Z = 50cm




Fig. 13 Lamp spectra FS40 and 0.127 mm, 0.191 mm, 0.254 mm cellulose acetate

Snowmass, Colorado




Fig. 14. Lamp spectra, Snowmass, Colorado

IRRADIANCE - $\text{mW} \cdot \text{M}^{-2} \cdot \text{nm}^{-1}$

100

9

8

7

6

5

4

3

2

1

0.1

0.09

0.08

0.07

0.06

0.05

0.04

0.03

0.02

0.01

0.009

0.008

0.007

0.006

0.005

0.004

0.003

0.002

0.001

0.0009

0.0008

0.0007

0.0006

0.0005

0.0004

0.0003

0.0002

0.0001

0.00009

0.00008

0.00007

0.00006

0.00005

0.00004

0.00003

0.00002

0.00001

0.000009

0.000008

0.000007

0.000006

0.000005

0.000004

0.000003

0.000002

0.000001

0.0000009

0.0000008

0.0000007

0.0000006

0.0000005

0.0000004

0.0000003

0.0000002

0.0000001

BZS 3-4 CLG 0377

20 W $Z = 50 \text{ cm}$

770523-35

$A_{21} = 6.9$

AGED 0h

VERTICAL
MOUNT

0h
3:30 PM

-16.5h

770524-15

$A_{21} = 4.3$

AGED 16.5h

$\lambda \text{ nm}$

Fig. 15. Lamp aging - BZS 20 CLG0377

NO 340 1510 DIETZEN (GRAPH PAPER)
F. CYLES X 10 DIV. (Hz) PER INCH

EUGENE




Fig. 16 Sylvania F15T8 CW lamp spectra - IRLSpec-S

EUGENE DIETZGEN CO
NEW YORK, N.Y.
5 CYCLES X 10 DIVISIONS PER INCH
SEMI-LOGARITHMIC




Fig. 17. Sylvania F15T8 CW lamp spectra IRLSpec-SO




Fig. 18. Sylvania F15T8 CW lamp spectra IRLSpec-D




Fig. 19. GE F15T8D lamp spectra IRLSpec-D




Fig. 20. Sun spectra - Beltsville, MD

NO 340 1510 DIETZGEN GRAPH PAPER
5 CYCLES X 1
EUGENE DIETZGEN CO




Fig. 21. Sun spectra - Beltsville, MD - IRLSpec-S




Fig. 22. Sun - Beltsville, MD - IRLSpec-D




Fig. 25. Regression IRLSpec-S vs. IRLMeter




Fig. 26. Rayonet F8T5 RPR 3000A⁰ spectra new cellulose acetate,
Z = 16.3 cm




Fig. 27. Rayonet F8T5 RPR 3000A° spectra cellulose acetate aged 6 hr.,
Z = 16.3 cm




Fig. 28. Rayonet F8T5 RPR 3000A⁰ spectra cellulose acetate aged 6 hr.
Z = 20 nm




Fig. 29. Rayonet F8T5 RPR 3000A^o and sun (Beltsville) spectra




Fig. 30. Spectra FS40 and Rayonet F8T5 RPR 3000A^o





Fig. 31. Spectra FS40 and cellulose acetate aged 0.5 hr. Z = 50 cm






Fig. 34. Spectra prescription eyeglasses and sunglasses

APPENDIX A

Figure

- A1a - A1e - Operate Program, IRLSpec-D
- A2 - Calibrate Program, IRLSpec-D
- A3a - A3b - Operate Program, IRLSpec-S
- A4a - A4d - Operate Program, IRLSpec-SO

Spectral Response Curves

- A5 - Standard Lamp Plotted Each nm
- A6 - FS40, FBZS40CLG1076, Sun, IRLSpec-S
- A7 - Filtered FS40 and FBZS40CLG1076, IRLSpec-S
- A8 - FS40, FBZS40CLG1076, Z = 73 cm, IRLSpec-D
- A9 - FBZS40 WLG, FBZS20 WLG, Vitalite, IRLSpec-D
- A10 - F15T8 CW Westinghouse, IRLSpec-D
- A11 - F15T8 CW Westinghouse, IRLSpec-SO
- A12 - F15T8 CW Sylvania, IRLSpec-D
- A13 - F15T8 GRO-LUX Sylvania F15T8 WW GE, IRLSpec-D
- A14 - F40R GE, F40IR Westinghouse, F15T8R GE, Z - 50 cm, IRLSpec-D
- A15 - F40IR Westinghouse, Z = 20 cm, IRLSpec-D
- A16 - F15T8BL Sylvania, F15T8BLB GE, IRLSpec-D
- A17 - F40 BL + FBZS40 CLG1076, IRLSpec-D
- A18 - F20T12/2021 Sample No. 3176-2, IRLSpec-D, IRLSpec-SO
- A19 - F20BL Phillips phosphor, F4T5BL conventional phosphor,
IRLSpec-S

AUTO START	0027 C	0075 .
A-IRP. CAL. SCAN	0028 A	0076
B-REV	0029 L	0077 D
C-CAL. WLGTH	0030 .	0078 A
D-SCAN SPECTRA	0031	0079 T
E-STOP	0032 W	0080 A
F-FORWARD	0033 L	0081 LINE
G-PRT. DATA	0034 G	0082 I
I-ACT. INT. 9	0035 T	0083 -
K-WLGTH. CHECK	0036 H	0084 A
L-STD. LAMP EQ	0037 LINE	0085 C
M-SUMS	0038 D	0086 T
N-CONT. READ	0039 -	0087 .
ENTER (R.MO.DAY#	0040 S	0088
2780120.00	0041 C	0089 I
RESET	0042 A	0090 N
	0043 H	0091 T
	0044	0092 .
FILE 0	0045 S	0093
TYPE 0	0046 P	0094 9
USED 746	0047 E	0095 LINE
MAX 1000	0048 C	0096 F
	0049 T	0097 -
	0050 R	0098 W
0000 PRNT*	0051 A	0099 L
0002 A	0052 LINE	0100 G
0003 -	0053 E	0101 T
0004 I	0054 -	0102 H
0005 R	0055 S	0103 .
0006 R	0056 T	0104
0007 .	0057 O	0105 C
0008	0058 P	0106 H
0009 C	0059 LINE	0107 E
0010 A	0060 F	0108 C
0011 L	0061 -	0109 F
0012 .	0062 F	0110 LINE
0013	0063 O	0111 L
0014 S	0064 R	0112 -
0015 C	0065 W	0113 S
0016 A	0066 A	0114 T
0017 H	0067 R	0115 D
0018 B	0068 D	0116 .
0019 -	0069 LINE	0117
0020 R	0070 G	0118 L
0021 E	0071 -	0119 A
0022 V	0072 P	0120 M
0023	0073 R	0121 F
0024 LINE	0074 T	0122
0025 C		0123 E
0026 -		0124 O

Fig. Ala. Operate Program IRLSpec-D (Cassette D4)

0125	LINE	0179	#REGS	0238	1
0126	M	0180	GOTO 0695	0239	ENTER↑
0127	-	0182	CLEAR	0240	0
0128	S	0183	6	0241	+*-
0129	U	0184	+*-	0242	LOAD
0130	M	0185	ENTER↑	0243	5
0131	S	0186	1	0244	WBYTE 4
0132	LINE	0187	FRMT 4	0246	READ 4
0133	N	0189	STOP	0248	ENTER↑
0134	-	0190	LBL	0249	4
0135	C	----	C	0250	+*-
0136	O	0192	LBL	0251	IF X=Y
0137	N	----	C	0252	GOTO 0256
0138	T	0194	PRNT*	0254	GOTO 0
0139	.	0196	N	0256	3
0140		0197	.	0257	WBYTE 4
0141	R	0198	L	0259	FOR A+F
0142	E	0199	.	0260	1
0143	A	0200		0261	STO B
0144	D	0201	C	0262	READ 4
0145	LINE	0202	A	0264	IF 0
0146	E	0203	L	0265	GOTO 0269
0147	N	0204	.	0267	GOTO 0
0148	T	0205	END*	0269	X=Y
0149	E	0206	5	0270	IF SFG 4
0150	R	0207	0	0271	GOTO 0530
0151		0208	STO F	0273	STO C
0152	Y	0209	SFG 2	0274	6
0153	R	0210	LBL	0275	WBYTE 4
0154	.	----	D	0277	FOR B+G
0155	M	0212	LBL	0278	READ 4
0156	O	----	D	0280	X=Y
0157	.	0214	FIX 2	0281	STO+ C
0158	D	0216	.	0282	NEXT B
0159	A	0217	0	0283	3
0160	Y	0218	1	0284	WBYTE 4
0161	#	0219	STO+ R000	0286	RCL C
0162	END*	0221	RCL R000	0287	RCL I A
0163	STOP	0223	PRINT	0289	-
0164	FIX 2	0224	FIX 4	0290	101X
0166	2	0226	1	0291	STO I A
0167	EEX	0227	STO A	0293	CFG 4
0168	6	0228	IF SFG 2	0294	RCL A
0169	+	0229	GOTO 0235	0295	2
0170	STO R000	0231	1	0296	4
0172	PRINT	0232	2	0297	3
0173	5	0233	0	0298	+
0174	MASK 4	0234	STO F	0299	ENTER↑
0176	1	0235	3	0300	1
0177	5	0236	9	0301	EEX
0178	1	0237	STO G	0302	3

Fig. Alb. Operate Program IRLSpec-D (Cassette D4) (cont'd)

0303	+		0365	END*		0418	RCL	I
0304	+		0366	2		0419	FIX	4
0305	PRINT		0367	8		0421	PRINT	
0306	IF CFG	8	0368	0		0422	ENTER↑	
0307	GOTO	0316	0369	STOP		0423	3	
0309	.		0370	FIX	0	0424	.	
0310	0		0372	PRINT		0425	0	
0311	1		0373	ENTER↑		0426	6	
0312	STO-	R000	0374	2		0427	+	
0314	GOTO	0317	0375	4		0428	FPNT*	
0316	NEXT	A	0376	9		0430	5	
0317	1		0377	-		0431	U	
0318	WBYTE	4	0378	STO	A	0432	N	
0320	CFG	8	0379	8		0433	5	
0321	RCL	A	0380	0		0434		
0322	8		0381	STO	F	0435	=	
0323	0		0382	CLEAR		0436	PRINT	
0324	-		0383	STO	I	0437	END*	
0325	IF +		0384	FOR	A+F	0438	GOTO	M
0326	GOTO	1	0385	RCL	A	0440	LBL	
0328	IF SFG	2	0386	2		----	A	
0329	GOTO	K	0387	4		0442	.	
0331	RETURN		0388	9		0443	0	
0332	LBL		0389	+		0444	1	
----	F		0390	2		0445	STO+	R000
0334	3		0391	2		0447	FIX	2
0335	WBYTE	4	0392	8		0449	RCL	R000
0337	RETURN		0393	.		0451	PRINT	
0338	LBL		0394	1		0452	FIX	3
----	E		0395	7		0454	1	
0340	CLEAR		0396	8		0455	STO	A
0341	WBYTE	4	0397	-		0456	1	
0343	RETURN		0398	9		0457	5	
0344	LBL		0399	Y↑		0458	0	
----	B		0400	STO	A	0459	STO	F
0346	1		0401	4		0460	4	
0347	WBYTE	4	0402	-		0461	0	
0349	1		0403	4		0462	STO	G
0350	RETURN		0404	Y↑		0463	1	
0351	LBL		0405	STO	J	0464	STO	E
----	I		0406	RCL	H	0465	5	
0353	FPNT*		0407	4		0466	WBYTE	4
0355	I		0408	-		0468	READ	4
0356	N		0409	+-		0470	ENTER↑	
0357	T		0410	0↑		0471	4	
0358	.		0411	RCL	J	0472	+-	
0359	9		0412	+		0473	IF X=Y	
0360			0413	RCL I	A	0474	GOTO	0478
0361	F		0415	+		0476	GOTO	0
0362	F		0416	ACC+		0478	3	
0363	0		0417	NEXT	A	0479	WBYTE	4
0364	M					0481	FOR	A+F

Fig. Alc. Operate Program IRLSpec-D (Cassette D4) (cont'd)

0482	READ	4	0546	RETURN	0603	END*
0484	IF 0		0547	LBL	0604	STOP
0485	GOTO	0489	----	L	0605	FIX 0
0487	GOTO	0	0549	LBL	0607	PRINT
0489	X=Y		----	L	0608	2
0490	IF SFG	4	0551	CLEAR	0609	4
0491	GOTO	0530	0552	1	0610	9
0493	CLEAR		0553	LD%GO	0611	-
0494	STO	C	0554	LBL	0612	STO A
0495	6		----	K	0613	PRNT*
0496	WBYTE	4	0556	RCL P047	0615	T
0498	FOR	B+G	0558	RCL R048	0616	0
0499	READ	4	0560	+	0617	
0501	X=Y		0561	LOG	0618	END*
0502	STO+	C	0562	ENTER↑	0619	STOP
0503	NEXT	B	0563	4	0620	PRINT
0504	3		0564	*	0621	2
0505	WBYTE	4	0565	+*-	0622	4
0507	1		0566	2	0623	9
0508	STO	B	0567	9	0624	-
0509	RCL	C	0568	6	0625	STO F
0510	STO+ 1	A	0569	.	0626	CLEAR
0512	+*-		0570	5	0627	STO I
0513	RCL	A	0571	+	0628	FOR A+F
0514	2		0572	FIX 2	0629	RCL A
0515	4		0574	PRNT*	0630	2
0516	9		0576	LINE	0631	4
0517	+		0577	L	0632	9
0518	ENTER↑		0578	I	0633	-
0519	EEX		0579	N	0634	RCL I A
0520	4		0580	E	0636	STO+ I
0521	*		0581		0637	NEXT A
0522	+		0582	A	0638	FIX 4
0523	PRINT		0583	T	0640	RCL I
0524	IF SFG 8		0584	PRINT	0641	PRNT*
0525	GOTO	0541	0585	END*	0643	M
0527	NEXT	A	0586	FIX 4	0644	I
0528	GOTO	0541	0588	CFG 2	0645	L
0530	PRNT*		0589	STOP	0646	L
0532	0		0590	LBL	0647	I
0533	V		----	M	0648	W
0534	E		0592	PRNT*	0649	A
0535	R		0594	S	0650	T
0536	L		0595	U	0651	T
0537	0		0596	M	0652	S
0538	A		0597		0653	.
0539	D		0598	F	0654	M
0540	END*		0599	R	0655	S
0541	1		0600	0	0656	0
0542	WBYTE	4	0601	M	0657	PRINT
0544	CFG	8	0602		0658	END*
0545	CFG	4				

Fig. Ald. Operate Program IRLSpec-D (Cassette D4) (cont'd)

0659	GOTO	M	0703	R	
0661	LBL		0704	E	
----	N		0705	S	
0663	CLEAR		0706	E	
0664	STO	C	0707	T	
0665	4		0708	END*	
0666	WBYTE	4	0709	GOTO	0182
0668	4		0711	LBL	
0669	0		----	G	
0670	STO	G	0713	LBL	
0671	1		----	G	
0672	STO	B	0715	RCL	P000
0673	FOR	B+G	0717	FIX	2
0674	READ	4	0719	PRINT	
0676	X*Y		0720	1	
0677	STO+	C	0721	STO	A
0678	NEXT	B	0722	1	
0679	RCL	C	0723	2	
0680	PAUSE		0724	0	
0681	PAUSE		0725	STO	F
0682	CLEAR		0726	FIX	4
0683	STO	C	0728	FOR	A+F
0684	GOTO	0671	0729	RCL I	A
0686	LBL		0731	RCL	A
----	0		0732	2	
0688	LBL		0733	4	
----	0		0734	9	
0690	.		0735	+	
0691	0		0736	ENTER†	
0692	1		0737	1	
0693	STO-	R000	0738	EEX	
0695	CLEAR		0739	3.	
0696	WBYTE	4	0740	+	
0698	1		0741	+	
0699	WBYTE	4	0742	PRINT	
0701	PRNT*		0743	NEXT	A
			0744	STOP	
			0745	END	

Fig. Ale. Operate Program IRLSpec-D (Cassette D4) (cont'd)

FILE	1	0051		0105	2	0161	+
TYPE	0	0052		0106	1	0162	RCL
USED	218	0053		0107	4	0163	LOG
MAX	800	0054	0	0108	STO	0164	-
		0055		0109	.	0165	STO I
		0056	A	0110	2	0167	PRINT
0000	PRNT	0057	H	0111	1	0168	NEXT
0002	C	0058	D	0112	4	0169	1
0003	A	0059		0113	6	0170	5
0004	L	0060	-	0114	STO	0171	0
0005	I	0061	1	0115	FOP	0172	ENTER↑
0006	B	0062	LINE	0116	RCL	0173	1
0007	R	0063	LINE	0117	2	0174	ENTER↑
0008	A	0064	LINE	0118	4	0175	0
0009	T	0065		0119	9	0176	+*-
0010	I	0066		0120	+	0177	RCDATA
0011	0	0067	H	0121	RCL	0178	1
0012	N	0068	0	0122	*	0179	5
0013	LINE	0069	.	0123	RCL	0180	0
0014	C	0070		0124	+	0181	ENTER↑
0015	0	0071	R	0125	10↑X	0182	1
0016	R	0072	U	0126	+*-	0183	ENTER↑
0017	R	0073	H	0127	RCL	0184	1
0018	E	0074	S	0128	+	0185	+*-
0019	C	0075	=	0129	10↑X	0186	RCDATA
0020	T	0076	?	0130	STO	0187	PRNT*
0021	I	0077	END*	0131	RCL	0189	C
0022	0	0078	STOP	0132	2	0190	A
0023	N	0079	STO	0133	4	0191	L
0024		0080	1	0134	9	0192	.
0025	W	0081	STO	0135	+	0193	
0026	I	0082	1	0136	RCL	0194	C
0027	L	0083	5	0137	-	0195	0
0028	L	0084	0	0138	.	0196	N
0029	LINE	0085	STO	0139	5	0197	S
0030	B	0086	.	0140	*	0198	T
0031	E	0087	0	0141	1↑X	0199	.
0032		0088	0	0142	STO-	0200	
0033	F	0089	3	0143	RCL	0201	I
0034	E	0090	5	0144	3	0202	N
0035	C	0091	+*-	0145	1	0203	LINE
0036	0	0092	STO	0146	-	0204	1
0037	R	0093	1	0147	RCL	0205	
0038	D	0094	.	0148	+	0206	T
0039	E	0095	3	0149	ENTER↑	0207	0
0040	D	0096	1	0150	*	0208	
0041		0097	7	0151	+*-	0209	1
0042	0	0098	STO	0152	10↑X	0210	5
0043	N	0099	2	0153	.	0211	0
0044	LINE	0100	.	0154	1	0212	END*
0045		0101	3	0155	+*-	0213	CLEAR
0046		0102	7	0156	*	0214	STO
0047	F	0103	6	0157	STO+	0215	CLEAR
0048	I	0104	STO	0158	RCL I	0216	LD%GO
0049	L			0160	RCL	0217	END
0050	E						

Fig. A2. Calibrate Program IRLSpec-D (Cassette D2)

DATE 780210
 E UW SPECTRA
 1780210.01

OPERATE PROGRAM IRLSpec-S

0000	CHA	0051	CHA	0101	I01	0151	RM	0201	04
0001	D	0052	U	0102	FLG	0152	11	0202	SP
0002	A	0053	W	0103	00	0153	-	0203	86
0003	T	0054		0104	IFER	0154	R	0204	71
0004	E	0055	S	0105	01	0155	IND	0205	R
0005	CHA	0056	P	0106	I06	0156	RM	0206	RM
0006	E	0057	E	0107	I01	0157	00	0207	41
0007	+	0058	C	0108	+	0158	=	0208	+
0008	1	0059	T	0109	RM	0159	10*	0209	L
0009	EXP	0060	R	0110	11	0160	R	0210	RM
0010	6	0061	A	0111	-	0161	IND	0211	42
0011	=	0062		0112	L	0162	SM	0212	=
0012	SM	0063	CHA	0113	IND	0163	01	0213	LOG
0013	10	0064	FIX5	0114	RM	0164	I02	0214	x
0014	E	0065	02	0115	00	0165	FIX5	0215	1
0015	SP	0066	COL	0116	=	0166	03	0216	=
0016	8a	0067	10	0117	10*	0167	COL	0217	90
0017	3	0068	CHA	0118	L	0168	08	0218	+
0018	7	0069		0119	IND	0169	+	0219	70
0019	0	0070	W	0120	SM	0170	01	0220	RM
0020	SM	0071	L	0121	01	0171	1	0221	13
0021	06	0072		0122	SM	0172	SM	0222	=
0022	FLG	0073		0123	03	0173	00	0223	71
0023	10	0074		0124	I02	0174	SM	0224	FIX5
0024	8	0075	I	0125	RM	0175	01	0225	02
0025	0	0076	R	0126	02	0176	SM	0226	COL
0026	SM	0077	R	0127	FIX5	0177	02	0227	07
0027	00	0078	A	0128	00	0178	RM	0228	L
0028	1	0079	D	0129	COL	0179	02	0229	RM
0029	4	0080	.	0130	03	0180	-	0230	63
0030	0	0081		0131	RM	0181	RM	0231	-
0031	SM	0082		0132	03	0182	06	0232	R
0032	01	0083	W	0133	FIX5	0183	=	0233	RM
0033	2	0084	L	0134	03	0184	IF-	0234	63
0034	5	0085		0135	COL	0185	00	0235	=
0035	0	0086		0136	08	0186	FLG	0236	LOG
0036	SM	0087		0137	+	0187	01	0237	x
0037	02	0088	I	0138	01	0188	I06	0238	2
0038	f9	0089	R	0139	1	0189	CHA	0239	=
0039	I06	0090	R	0140	SM	0190	0	0240	90
0040	CHA	0091	A	0141	02	0191	CHA	0241	+
0041	2	0092	D	0142	RM	0192	I02	0242	70
0042	CHA	0093	.	0143	02	0193	LF	0243	RM
0043	I02	0094	CHA	0144	FIX5	0194	RM	0244	14
0044	.	0095	+	0145	00	0195	02	0245	=
0045	0	0096	01	0146	COL	0196	-	0246	FIX5
0046	1	0097	I06	0147	03	0197	RM	0247	02
0047	SM	0098	CHA	0148	I06	0198	07	0248	COL
0048	10	0099	7	0149	I01	0199	=	0249	07
0049	RM	0100	CHA	0150	+	0200	IF+	0250	LF
0050	10								

Fig. A3a

OPERATE PROGRAM IRLSpec-S

(cont'd)

0251 E	0301 FLG	0352 CHA	0401 IND	0451
0252 FLG	0302 02	0353 +	0402 RM	0452 F
0253 04	0303 L	0354 RM	0403 03	0453 R
0254 SP	0304 IND	0355 15	0404 SM	0454 0
0255 00	0305 RM	0356 =	0405 02	0455 M
0256 CM	0306 00	0357 FIX5	0406 R	0456 CHA
0257 02	0307 X	0358 05	0407 IND	0457 LF
0258 CHA	0308 L	0359 COL	0408 RM	0458 E
0259 S	0309 IND	0360 16	0409 03	0459 SM
0260 T	0310 RM	0361 LF	0410 SM	0460 08
0261 A	0311 01	0362 SP	0411 02	0461 CHA
0262 R	0312 =	0363 0d	0412 1	0462 T
0263 T	0313 SM	0364 FLG	0413 SM	0463 0
0264	0314 02	0365 07	0414 03	0464 CHA
0265 I	0315 R	0366 CHA	0415 RM	0465 E
0266 N	0316 IND	0367 S	0416 03	0466 LF
0267 T	0317 RM	0368 U	0417 -	0467 SM
0268 .	0318 00	0369 M	0418 1	0468 09
0269 9	0319 X	0370	0419 5	0469 FLG
0270	0320 R	0371 F	0420 =	0470 11
0271 A	0321 IND	0372 R	0421 X	0471 L
0272 T	0322 RM	0373 0	0422 2	0472 IND
0273 ?	0323 01	0374 M	0423 -	0473 RM
0274 CHA	0324 =	0375 CHA	0424 RM	0474 08
0275 E	0325 SM	0376 E	0425 04	0475 FIX5
0276 IFE	0326 02	0377 SM	0426 =	0476 03
0277 12	0327 1	0378 03	0427 IF-	0477 COL
0278 2	0328 SM	0379 CHA	0428 03	0478 12
0279 8	0329 00	0380 T	0429 LF	0479 R
0280 0	0330 SM	0381 0	0430 RM	0480 IND
0281 0	0331 01	0382 CHA	0431 02	0481 RM
0282 FLG	0332 RM	0383 E	0432 FIX5	0482 08
0283 12	0333 01	0384 SM	0433 03	0483 FIX5
0284 -	0334 -	0385 04	0434 0	0484 03
0285 1	0335 7	0386 CM	0435 GT	0485 COL
0286 7	0336 9	0387 02	0436 07	0486 10
0287 0	0337 =	0388 RM	0437 SP	0487 LF
0288 =	0338 IF-	0389 03	0438 06	0488 1
0289 +	0339 02	0390 +	0439 3	0489 SM
0290 2	0340 RM	0391 2	0440 0	0490 08
0291 =	0341 02	0392 +	0441 0	0491 RM
0292 SM	0342 FIX5	0393 1	0442 SM	0492 08
0293 01	0343 03	0394 5	0443 06	0493 -
0294 +	0344 COL	0395 =	0444 GT	0494 RM
0295 1	0345 08	0396 SM	0445 10	0495 09
0296 0	0346 CHA	0397 03	0446 CHA	0496 =
0297 0	0347	0398 FLG	0447 M	0497 IF-
0298 =	0348 S	0399 03	0448 E	0498 11
0299 SM	0349 U	0400 L	0449 M	0499 EP
0300 00	0350 N		0450 .	0500 00
	0351 S			

Fig. A3b

AUTO START	0023 C	0074 LINE	0123 LINE
A SPECTRA	0024	0075 J	0124 O
B REVERSE	0025 S	0076	0125
C S L CAL.	0026	0077 R	0126 O
E STOP	0027 L	0078 U	0127 P
F FORWARD	0028	0079 N	0128 E
G PRINT DATA	0029 C	0080	0129 R
I INTEGRAL	0030 A	0081 W	0130 A
J RUN WLGTH	0031 L	0082 L	0131 T
K WLGTH EQUAT.	0032 .	0083 G	0132 E
L STD LAMP EQUA	0033 LINE	0084 T	0133
M SUM	0034 E	0085 H	0134 LINE
O OPERATE	0035	0086 LINE	0135 E
ENTER YR.MO.DAY#	0036 S	0087 K	0136 H
	0037 T	0088	0137 T
	0038 O	0089 W	0138 E
3780120.00	0039 P	0090 L	0139 R
READY	0040 LINE	0091 G	0140
HIT A C OR J	0041 F	0092 T	0141 Y
	0042	0093 H	0142 R
FILE 0	0043 F	0094	0143 .
TYPE 0	0044 O	0095 E	0144 M
USED 720	0045 R	0096 O	0145 O
MAX 800	0046 W	0097 U	0146 .
	0047 A	0098 A	0147 D
0000 PRNT*	0048 R	0099 T	0148 A
0002 A	0049 D	0100 .	0149 Y
0003	0050 LINE	0101 LINE	0150 #
0004 S	0051 G	0102 L	0151 LINE
0005 P	0052	0103	0152 END*
0006 E	0053 P	0104 S	0153 STOP
0007 C	0054 R	0105 T	0154 3
0008 T	0055 I	0106 D	0155 EEX
0009 R	0056 H	0107	0156 6
0010 A	0057 T	0108 L	0157 +
0011 LINE	0058	0109 A	0158 FIX 2
0012 B	0059 D	0110 M	0160 PRINT
0013	0060 A	0111 P	0161 STO P000
0014 R	0061 T	0112	0163 2
0015 E	0062 A	0113 E	0164 5
0016 V	0063 LINE	0114 O	0165 4
0017 E	0064 I	0115 U	0166 0
0018 R	0065	0116 A	0167 STO B
0019 S	0066 I	0117 LINE	0168 2
0020 E	0067 N	0118 M	0169 5
0021	0068 T	0119	0170 1
0022 LINE	0069 E	0120 S	0171 0
	0070 G	0121 U	0172 STO D
	0071 R	0122 M	0173 FIX 4
	0072 A		
	0073 L		

Fig. A4a. Operate Program IRLSpec-SO (Cassette B)

0175	1		0229	STO+	R000	0290	LBL		0342	9
0176	5		0231	RCL	R000	----	I		0343	Y↑X
0177	1		0233	PRINT		0292	LBL		0344	STO
0178	#REGS		0234	FIX	4	----	I		0345	4
0179	LBL		0236	3		0294	PRNT%		0346	+
----	0		0237	WBYTE	4	0296	I		0347	4
0181	LBL		0239	1		0297	N		0348	Y↑X
----	0		0240	STO	A	0298	T		0349	STO
0183	CFG	8	0241	IF SFG	4	0299	.		0350	RCL
0184	CFG	3	0242	GOTO	0248	0300			0351	4
0185	PRNT%		0244	1		0301	F		0352	-
0187	R		0245	2		0302	R		0353	+*-
0188	E		0246	0		0303	0		0354	e↑X
0189	A		0247	STO	F	0304	M		0355	RCL
0190	D		0248	FOR	A+F	0305			0356	*
0191	Y		0249	START	4	0306	?		0357	RCL I
0192	LINE		0251	READ	4	0307	END%		0359	+
0193	H		0253	1		0308	2		0360	ACC+
0194	I		0254	EEX		0309	8		0361	NEXT
0195	T		0255	3		0310	0		0362	RCL
0196			0256	+		0311	STOP		0363	PRNT%
0197	A		0257	STO- I	A	0312	FIX	0	0365	LINE
0198			0259	RCL I	A	0314	PRINT		0366	I
0199	C		0261	+*-		0315	FIX	4	0367	N
0200			0262	10↑X		0317	ENTER↑		0368	T
0201	0		0263	STO I	A	0318	2		0369	E
0202	R		0265	RCL	E	0319	4		0370	G
0203			0266	-		0320	9		0371	R
0204	J		0267	IF -		0321	-		0372	A
0205	LINE		0268	CLEAR		0322	STO	A	0373	L
0206	END%		0269	RCL	A	0323	8		0374	
0207	STOP		0270	2		0324	0		0375	9
0208	LBL		0271	4		0325	STO	F	0376	
----	A		0272	9		0326	CLEAR		0377	=
0210	LBL		0273	+		0327	STO	I	0378	PRINT
----	A		0274	ENTER↑		0328	FOR	A→F	0379	LINE
0212	1		0275	1		0329	RCL	A	0380	LINE
0213	ENTER↑		0276	EEX		0330	2		0381	END%
0214	HNDSK	4	0277	3		0331	4		0382	ENTER↑
0216	2		0278	*		0332	9		0383	3
0217	WBYTE	4	0279	+		0333	+		0384	.
0219	1		0280	PRINT		0334	2		0385	0
0220	ENTER↑		0281	IF SFG	8	0335	2		0386	6
0221	0		0282	GOTO	B	0336	8		0387	+
0222	+*-		0284	NEXT	A	0337	.		0388	PRNT%
0223	LOAD		0285	IF SFG	4	0338	1		0390	S
0224	FIX	2	0286	GOTO	K	0339	7		0391	U
0226	.		0288	GOSUB	B	0340	8		0392	N
0227	0					0341	+		0393	S
0228	1									

Fig. A4b. Operate Program IRLSpec-SO (Cassette B) (cont'd)

0394		0450	GOSUB	B	0511	PAUSE	
0395	=	0452	GOTO	0	0512	PAUSE	
0396		0454	LBL		0513	FOR	A+F
0397	PRINT	----	J		0514	READ	4
0398	LINE	0456	LBL		0516	1	
0399	LINE	----	J		0517	EEP	
0400	END*	0458	5		0518	3	
0401	GOTO	0459	0		0519	-	
0403	LBL	0460	STO	F	0520	STO+ I	A
----	K	0461	SFG	4	0522	+*-	
0405	LBL	0462	GOTO	A	0523	RCL	A
----	K	0464	LBL		0524	2	
0407	RCL	----	F		0525	4	
	R004	0466	2		0526	9	
0409	RCL				0527	+	
	R005	0467	WBYTE	4	0528	ENTER↑	
0411	-	0469	3		0529	1	
0412	LOG	0470	WBYTE	4	0530	EEP	
0413	+*-	0472	RETURN		0531	3	
0414	ENTER↑	0473	LBL		0532	+	
0415	3	----	E		0533	+	
0416	.	0475	2		0534	FIX	4
0417	4	0476	WBYTE	4	0536	PRINT	
0418	+	0478	CLEAR		0537	IF SFG	6
0419	2	0479	WBYTE	4	0538	GOTO	B
0420	5	0481	RETURN		0540	NEXT	A
0421	3	0482	LBL		0541	LBL	
0422	.	----	C		----	B	
0423	5	0484	LBL		0543	LBL	
0424	+	----	C		----	B	
0425	FIX	0486	.		0545	1	
	2	0487	0		0546	ENTER↑	
0427	PRINT	0488	1		0547	HNDST	4
0428	RCL	0489	STO+	R000	0549	2	
	R047	0491	RCL	R000	0550	WBYTE	4
0430	RCL	0493	FIX	2	0552	PAUSE	
	R048	0495	PRINT		0553	0	
0432	-	0496	1		0554	WBYTE	4
0433	LOG	0497	STO	A	0556	PAUSE	
0434	+*-	0498	1		0557	1	
0435	ENTER↑	0499	5		0558	WBYTE	4
0436	3	0500	0		0560	PAUSE	
0437	.	0501	STO	F	0561	5	
0438	5	0503	1		0562	WBYTE	4
0439	*	0504	ENTER↑		0564	READ	4
0440	2	0505	HNDST	4	0566	READ	4
0441	9	0507	3		0568	RCL	B
0442	6	0508	WBYTE	4	0569	+	
0443	.	0510	PAUSE		0570	IF -	
0444	5						
0445	+						
0446	PRINT						
0447	FIX						
	4						
0449	CFG						
	4						

Fig. A4c. Operate Program IRLSpec-SO (Cassette B) (cont'd)

0571 GOTO	0566	0626 2	0676 LBL
0573 1		0627 2	---
0576 READ	4	0628	0678 LBL
0578 RCL	D	0629 -	---- G
0579 +		0630 STO	0680 1
0580 IF -		0631 PRNT*	0681 STO
0581 GOTO	0576	0633 T	0682 1
0583 0		0634 0	0683 5
0584 WBYTE	4	0635 END*	0684 0
0586 IF SFG	3	0636 STOP	0685 STO
0587 STOP		0637 PRINT	0686 FIX
0588 IF SFG	8	0638 2	0688 FOR
0589 GOTO	0712	0639 4	0689 RCL I
0591 RETURN		0640 9	0691 RCL
0592 LBL		0641 -	0692 2
---- M		0642 STO	0693 4
0594 LBL		0643 CLEAR	0694 9
---- M		0644 STO	0695 +
0596 PRNT*		0645 FOR	0696 ENTER
0598 R		0646 RCL I	0697 1
0599 E		0648 STO+	0698 EEX
0600 C		0649 NEXT	0699
0601 0		0650 FIX	0700
0602 R		0652 RCL	0701 +
0603 D		0653 PRNT*	0702 PRINT
0604 ?		0655 M	0703 NEXT
0605		0656 I	0704 STOP
0606 0		0657 L	0705 LBL
0607 R		0658 L	---- L
0608 LINE		0659 I	0707 LBL
0609 S		0660 W	---- L
0610 U		0661 A	0709 CLEAR
0611 M		0662 T	0710 1
0612		0663 T	0711 LD&GO
0613 F		0664 S	0712 .
0614 R		0665 /	0713 0
0615 0		0666 M	0714 1
0616 M		0667 S	0715 STO-
0617 END*		0668 0	0717 GOTO
0618 CLEAR		0669	0719 END
0619 STOP		0670 PRINT	
0620 IF 0		0671 LINE	
0621 GOTO	0	0672 LINE	
0623 FIX	0	0673 END*	
0625 PRINT		0674 GOTO	M

Fig. A4d. Operate Program IRLSpec-SO (Cassette B) (cont'd)




Fig. A5 - Standard Lamp Plotted Each nm

SPECTRAL IRRADIANCE - $\text{mW} \cdot \text{m}^{-2} \cdot \text{nm}^{-1}$




Fig. A6 - FS40, FBZS40CLG1076, Sun, IRLSpec-S




Fig. A7 - Filtered FS40 and FBZS40CLG1076, IRLSpec-S

Z 73 CM

ENERGY
270-340
521.4
1,883.4
1,109.
853.3

7-29-77

7-29-77

SPECTRAL IRRADIANCE - $\text{mW} \cdot \text{m}^{-2} \cdot \text{nm}^{-1}$




Fig. A8 - FS40, FBZS40CLG1076, Z = 73 cm, IRLSpec-D

NO. 440-1000 DITZGEN 1000-1000000
 1000-1000000 1000-1000000
 1000-1000000 1000-1000000




Fig. A10 - F15T8 CW Westinghouse, IRLSpec-D

Q. 11

ND 340 1510 DIEZIGEN ...




Fig. All - F15T8 CW Westinghouse, IRLSpec-60

NO. 340 1510 DIETZEN GRAPH PAPER
 100 MILLIGRAM PER INCH
 5 CYCLES X 10 DIVISIONS PER INCH




Fig. A12 - F15T8 CW Sylvania, IRLSpec-D

SPECTRAL IRRADIANCE - $\text{mW} \cdot \text{m}^{-2} \cdot \text{nm}^{-1}$




Fig. A13 - F15T8 GRO-LUX Sylvania F15T8 WW GE, IRLSpec-D

IRLSpec-D UV-B (2nm BW)
IRL-AMRI-ARS-USDA




Fig. A14 - F40R GE, F40IR Westinghouse, F15T8R GE, $Z = 50 \text{ cm}$,
IRLSpec-D




Fig. A15 - F40IR Westinghouse, Z = 20 cm, IRLSpec-D

NO. 340 11510 DIETZGEN GRAPH PAPER
SEMI-LOGARITHMIC
5 CYCLES X 10 DIVISIONS PER INCH




Fig. A16 - F15T8BL Sylvania, F15T8BLB GE, IRLSpec-D

SPECTRAL IRRADIANCE - $\text{mW} \cdot \text{m}^{-2} \cdot \text{nm}^{-1}$




Fig. A17 - F40 BL + FBZS40 CLG1076, IRLSpec-D

SPECTRAL IRRADIANCE - $\text{mW} \cdot \text{m}^{-2} \cdot \text{nm}^{-1}$




Fig. A18 - F20T12/2021 Sample No. 3176-2, IRLSpec-D, IRLSpec-SO

EUGENE DIETZGEN CO.
MADE IN U.S.A.

NO. 340-LS10 DIETZGEN GRAPH PAPER
SEMI-LOGARITHMIC
5 CYCLES X 10 DIVISIONS PER INCH

SPECTRAL IRRADIANCE - $\text{mW} \cdot \text{m}^{-2} \cdot \text{nm}^{-1}$




Fig. A19 - F20BL Phillips phosphor, F4T5BL conventional phosphor, IRLSpec-S

NATIONAL AGRICULTURAL LIBRARY



1022397545